

# Time Spectrometer APV85G4

MADE IN JAPAN

5GHz 4CH 10bit ADC Digital waveform processing

VME / UNIT

This is a spectrometer for time analysis that carried high-speed 5GHz ADC every CH. Unified functions of several radiation measurement modules which is necessary for time analysis such as Differential CFD, Delay TAC, MCA etc. The pre-amp signal from detector is sampled to wave form in high-speed ADC and time analyzed it in FPGA and transfer the operation result to a PC by Ethernet.

High time resolution  
High throughput




## Features

- ◆ **ADC** 4CH, 5GHz Sampling, 10bit resolution
- ◆ **Time Resolution** Course: 200ps | Fine: 0.78ps, LSB
- ◆ **Throughput** 1Mcps or more /CH
- ◆ **Interface** TCP/IP(Gigabit-Ether) Data transfer 20MByte/sec or more
- ◆ **Analyze Mode** Differential-TDC, List(TDC+QDC etc.), Waveform, Histogram etc.
- ◆ **Function** Digital CFD(WALK, THRESHOLD, LLD, ULD), PSA, coincidence
- ◆ **Usage Example** Positron lifetime measurement,  
Discrimination of  $\gamma$  and neutron by liquid scintillator.  
High-speed and high resolution for scintillator such as  $\text{LaBr}_3(\text{Ce})$  or  $\text{CeBr}_3$

### Example ①

#### [Positron lifetime measurement]




Example: Time resolution: Under 160ps  
( $\text{BaF}_2$  scintillator vs  $\text{BaF}_2$  scintillator), Sample: Silicon

Combination  
measurement of  
different modules


Unified functions  
of each modules

### Example ②


#### [Discrimination of $\gamma$ and neutron by scintillator]



## Waveform sampling



Pulsar



$\text{BaF}_2$  @ 511 keV ( $^{22}\text{Na}$ )

\*Images is for illustration purpose.  
\*Please note that contents may change without prior notice.

## Specifications

Analog signal input type	PMT anode signal, Fast-NIM signal etc.
Analog signal input range	$\pm 0.5\text{V}$ (Input impedance : 50 $\Omega$ ) (SMA connector x 4) *Customizable maximum $\pm 4\text{V}$ capable
Analog Offset adjustment	$\pm 0.5\text{V}$ (Customizable 20mV to $\pm 4\text{V}$ capable)
Analog signal rise time	Under 0.5ns
Analog gain switch	x 1 / x 3 *Customizable ATT capable
Time differential range	12.8ps to 3.3 $\mu\text{sec}$ *0.78ps/ch to 200ps/ch x 16384 ch
Outside In/output signal terminal	LEMO connector x 6 (TTL Level)
Electricity consumption	+5V (9A), +12V (1A), -12V(1A)
Environmental conditions	Temperature 5 to 25 $^{\circ}\text{C}$
Dimension and Weight	VME1width 20mm(W) x 262mm(H) x 187mm(D), 520g
Accessory	Application and Manual

### List Data Example (1 event: 80bit)

80 15 13 0  
TDC[63..0] CH#[1..0] QDC[13..0]

# TechnoAP

Design and fabrication of electronic circuit associated  
with measurement control and radiation measurement

TechnoAP Co., Ltd.

+81-29-350-8011

+81-29-352-9013

2976-15 Mawatari, Hitachinaka-shi, Ibaraki, 312-0012, Japan

<http://www.techno-ap.com>

[order@techno-ap.com](mailto:order@techno-ap.com)

Updated on 2017/03/29