

Digital Pulse Processor APV8108-14

Instruction Manual

Version 1.4
September 2018

TechnoAP Co., Ltd.

ADD : 2976-15 Mawatari, Hitachinaka-shi, Ibaraki, Japan

ZIP Code: 312-0012

TEL: 029-350-8011

FAX: 029-352-9013

URL: <http://www.techno-ap.com>

e-mail: order@techno-ap.com

- Contents -

1.	Safety Precautions / Disclaimer.....	3
2.	Overview	4
2. 1.	Overview.....	4
2. 2.	Specifications	5
2. 3.	Revision History	5
3.	Appearance	6
3. 1.	Appearance.....	6
4.	Setup	7
4. 1.	Instration of the application.....	7
4. 2.	Connection	7
4. 3.	Setup of the network.....	8
5.	Screen of the application.....	10
5. 1.	Startup screen	10
5. 2.	Config tab.....	13
5. 3.	file tab.....	22
5. 4.	wave tab.....	24
5. 5.	spectrum tab.....	26
5. 6.	timespectrum tab	28
6.	Measurement.....	29
6. 1.	Energy Spectrum measurement	29
6. 2.	List measurement	33
6. 3.	Time Spectrum measurement.....	36
7.	File	41
7. 1.	Histogram data file	41
7. 2.	Wave data file.....	43
7. 3.	List data file	45
8.	Command	46
8. 1	Overview	46
8. 2	Format of command	47
8. 3	Type of command.....	48
8. 4	List of command.....	51
8. 5	Explanation of command.....	54
8. 6	Setting command at startup and config.....	64
9.	End.....	77

1. Safety Precautions / Disclaimer

Thank you very much for purchasing the digitizer APV8108-14 (hereinafter "This board") of TechnoAP Co., Ltd. (hereinafter "We"). Please read this "Safety Precautions / Disclaimer" before using this device, be sure to observe the contents, and use it correctly.

We are not responsible for any damage caused by abnormality of device, detector, connected device, application, damage to failure, other secondary damage, even if accident caused by using this device.

Prohibited matter

- This device cannot be used for applications requiring special quality and reliability related to human life, accident.
- This device cannot be used in places with high temperature, high humidity and high vibration.
- Do not apply a power supply that exceeds the rating.
- Do not turn the power on while other metals are in contact with the board surface.

Note

- If there is smoking or abnormal heat generation in this device, turn off the power immediately.
- This board may not work properly in noisy environments.
- Be careful with static electricity.
- The specifications of this board and the contents of the related documents are subject to change without notice.

Warranty policy

The warranty conditions of "our product" are as follows.

Warranty period	One year from date of purchase.
Guarantee contents	Repair or replacement will be carried out in case of breakdown even though you have used correctly according to this instruction manual within the warranty period
Out of warranty	We do not warranty if the cause of the failure falls under any of the following. 1. Failure or damage due to misuse or improper repair or modification or disassembly. 2. Failure and damage due to falling etc. 3. Breakdown / damage in harsh environments (high temperature / high humidity, under zero, condensation etc.). 4. Causes other than the above, other than "our products". 5. Consumables.

2. Overview

2. 1. Overview

The APV8108-14 is a waveform analysis board that uses a high-speed, high-resolution ADC. Real-time analysis at 1 GHz with FPGA, high speed processing without dead time by signal processing is realized with high time resolution and high throughput. All ADCs operate synchronously with a 1 GHz clock and can be used for signal analysis from multiple high-speed scintillation detectors. In addition, it supports synchronous processing of multiple boards, and can easily be extended to multi-channel analysis.




Fig. 1 Function

This manual describes the software for measuring and controlling this device.

* "list" and "event" in the text are equivalent.

* "hist" and "spectrum" in the text are equivalent.

2. 2. Specifications


- (1) Analog Input
- Number of channel : 8 channel
 - Input range : ± 1 V
 - Input impedance : 50 Ω
- (2) ADC
- Sampling frequency : 1 GHz
 - Resolution : 14-bit
 - SNR : 68.3 dBFS @ 605 MHz
- (3) Performance
- QDC output : 2 Mcps and more
 - Time resolution : 3.90625 ps
- (4) MCA
- Measurement mode : Wave mode, Histogram mode, List mode, List-Common mode
 - Transfer rate : Appx. 20 Mbyte per second
- (5) Interface
- LAN : Ethernet TCP/IP 1000 Base-T (at List data acquisition)
UDP (When sending / receiving config data, receiving status data)
- (6) Form
- VME type : APV8108-14
- (7) Consumption current
- + 5 V : 6.0 A (Max.)
 - + 12 V : 1.0 A (Max.)
- (8) Application
- OS : Windows 7 version or later, 32-bit or 64-bit
 - Screen resolution : Recommend HD (1366 x 768) or more

2. 3. Revision History

Version 1.0	May 2018	First edition
Version 1.1	June 2018	Correction of errors
Version 1.2	July 2018	Correction of errors
Version 1.3	August 2018	Correction of errors
Version 1.4	September 2018	Correction of errors

3. Appearance

3. 1. Appearance



Pic. 1: APV8108-14

1	LED	P: Power ON, V: Not used. E: Not used.
2	CH1 to CH8	LEMO connector for signal input. Input range: ± 1 V, input impedance: 50 Ω .
3	SYNC-O	LEMO connector for synchronous timing signal output. Outputs a timing signal to adjust the time between boards.
4	SYNG-I	LEMO connector for synchronous timing signal input. Input a timing signal to adjust the time between boards.
5	CLK-O	LEMO connector for external clock signal output. Outputs a 25 MHz TTL signal.
6	CLK-I	LEMO connector for external clock signal input. It can be operated using an external clock. Turn on the power after inputting the 25 MHz TTL signal.
7	VETO	LEMO connector for external signal input. Disable data acquisition during "High".
8	GATE	LEMO connector for external gate signal input. Input TTL signal. Data acquisition is enabled while the input is "High".
9	LAN	RJ45 connector for Ethernet cables. 1000Base-T.

NOTE: Connect the SYNC-O and SYNC-I by using a cable.

4. Setup

4. 1. Instration of the application

The application for APV8108-14 (hereinafter referred to as “the application”) runs on Windows. When using it, it is necessary to install this application's EXE (executable) file and National Instruments' LabVIEW Run-Time Engine on the PC used for measurement.

Installation of this application is performed by the installer included in the attached CD. The installer includes an EXE (executable) file and the LabVIEW Run-Time Engine, which can be installed at the same time.


The installation procedure is as follows.

1	Log in to Windows with administrator privileges.
2	"Setup.exe" in the "Installer" folder in the supplied CD-ROM. Proceed with the installation interactively. The default installation destination is "C: ¥ TechnoAP".
3	"Start button"-"TechnoAP"-"APV8108-TOTAL-FALL" is executed.

To uninstall, select "APV8108-TOTAL-FALL" from "Add or Remove Programs" to remove.


4. 2. Connection

1	Connect this device and PC with an Ethernet cable. Use a cross cable depending on the PC. Use a switching hub when using a hub.
2	SYNC-O terminal-SYNC-I terminal connection Connect the APV8108-14 SYNC-O and SYNC-I pins together.



Pic. 2: SYNC-I / O terminal connection example

4. 3. Setup of the network

1	Turn on the power of the PC and change the network information of the PC.	
	IP Address	192.168.10.2 *Any value except 192.168.10.128
	Sub-net mask	255.255.255.0
	Default gateway	192.168.10.1
2	Power on the VME rack. Do not operate anything for 10 seconds after power on.	
3	Check the communication connection between the PC and this device. Execute the ping command at the Windows command prompt to check whether the device and PC can be connected. The IP address of this device is on the board. The network information of this device at the time of shipment from the factory is as follows.	
	IP Address	192.168.10.125 *The following is an example at 192.168.10.128.
	Sub-net mask	255.255.255.0
	Default gateway	192.168.10.1
	<p>> ping 192.168.10.128</p>  <pre> C:\WINDOWS\system32\cmd.exe Microsoft Windows [Version 10.0.14393] (c) 2016 Microsoft Corporation. All rights reserved. C:\Users\Administrator>ping 192.168.10.128 192.168.10.128 に ping を送信しています 32 バイトのデータ: 192.168.10.128 からの応答: バイト数 =32 時間 <1ms TTL=128 192.168.10.128 からの応答: バイト数 =32 時間 <1ms TTL=128 192.168.10.128 からの応答: バイト数 =32 時間 <1ms TTL=128 192.168.10.128 からの応答: バイト数 =32 時間 <1ms TTL=128 192.168.10.128 の ping 統計: パケット数: 送信 = 4、受信 = 4、損失 = 0 (0% の損失)、 ラウンド トリップの概算時間 (ミリ秒): 最小 = 0ms、最大 = 0ms、平均 = 0ms C:\Users\Administrator> </pre>	
Pic. 3: Communication connection check, ping command execution		
4	Launch this application on a PC.	
	*When starting this application, an error message of the content that failed to connect with the device may be displayed. The main causes are as follows.	
	1. The port definition in the "System" section in the configuration file "config.ini" has an incorrect value. In particular, "DevConfigPort = 4660", "DevDataPort = 24", "SubnetMask =" 255.255.255.0 "", "Gateway =" 192.168.10.1 "", "ChNumber = 16" are important.	
	2. The LAN cable on the PC side is not fully inserted.	
	3. The LAN cable on the device side is not fully inserted.	
	4. The power of this device remains OFF or the LAN cable is broken.	
	5. The network setting on the PC side is DHCP.	
6. The network settings on the PC side are not configured with private addresses (192.168.10.2 to 255 excluding 192.168.10.128).		
7. The power saving mode of the PC is functioning.		

8. Wireless LAN of PC is enabled.


If the above cause does not start up properly, please try the following method.

After confirming the cable connection etc., restart this application.

5. Screen of the application

5. 1. Startup screen

The following opening screen is displayed when you launch "Start"- "TechnoAP"- "APV8108-14"



Pic. 4: Startup screen

The contents of each item are as follows.

Menu

It consists of "File", "Edit", "Calibration", "Config", "Clear", "Start" and "Stop".

File - open config	Load configuration file
File - save config	Save current settings to a file
File - save histogram	Save current histogram data to file
File - save wave	Save current waveform data to file
File - save image	Save this application screen as PNG format image
File - quit	Finish
Edit - copy setting of CH1	Reflects the setting of CH1 in the "CH" tab to the settings of all other CHs
Edit - IP configuration	Change IP address of display device
Calibration	Run calibration. Execute if there is a disturbance in the wave waveform
Help	Open the description of the current tab in pdf format
Config	Send all settings to this device
Clear	Initialize the histogram data in this device
Start	Send measurement start to this device
Stop	Send measurement stop to this device

Tab

config	Settings related to this device settings and measurement
file	Settings for saving waveform and list data
save	Display of input waveform, CFD waveform, QDC waveform
spectrum	Histogram display
timespectrum	Display time difference spectrum from time information of list data

CH section

Display the status of each CH.

output count	The total number of output events
output rate (cps)	Number of output events per second
deadtime (%)	Ratio of dead time

ROI section

Display the calculated result between ROIs.


peak (ch)	ch of the maximum count
centroid (ch)	Center value (ch) calculated from the sum of all counts
peak (count)	Maximum count
gross (count)	Sum of counts between ROIs
gross (cps)	CPS of counts between ROIs
net (count)	Sum of counts minus background between ROIs
net (cps)	CPS of the count minus the background between ROIs
FWHM (ch)	Half width (ch)
FWHM (%)	Half width (%). Half-width ÷ ROI defined energy × 100
FWHM	Half width
FWTM	1/10 width

Other

device	Select the device to be measured
IP address	It is defined in the configuration file, and the IP address of the device selected by "Module" is displayed.
memo	You can enter a memo for image storage.
acq. LED	Blink while measurement
save LED	Blink while saving list data
error LED	Lights up when an error occurs
mode	Display "hist", "wave", "list" or "list-com". Please note that the above modes may not be available depending on the configuration of the option.
measurement mode	Display "real time" and "live time"
measurement time	Set measurement time

real time	Real time of the effective first CH (actual measurement time). It is equal to measurement time at the end of measurement
live time	Effective first CH live time (effective measurement time). The value of real time minus the value of dead time.
file size (Byte)	Display the file capacity (Byte) while saving event data.
sampling	Display sampling frequency of display device

5. 2. Config tab




Pic. 5-1: config tab, upper part

The contents of each item are as follows.


CH enable	<p>CH availability.</p> <p>Normally, set all channels to enable (pressed).</p>
signal type	<p>Select the type of input waveform.</p> <p>Set to "fast sig" when inputting NIM signals or timing signals.</p> <p>For others, set "nomal sig".</p>
signal delay	<p>The input signal is delayed internally in this device.</p> <p>The maximum delay time is 2 μs.</p>
polarity	<p>Select the polarity of the input signal from "pos" for positive polarity and "neg" for negative polarity.</p>
baseline restorer filter	<p>Set the time constant of baseline restorer.</p> <p>Set from Ext (without AutoBLR), Fast, 4 μs, 85 μs, 129 μs, and 260 μs.</p> <p>Normally set to 85 μs.</p>
threshold	<p>Set the threshold for waveform acquisition of the input signal. The unit is digit. The setting range is 0 to 8191. While watching the "raw" waveform in wave mode, set the value higher than the noise level.</p> <div style="text-align: center;"> <p>threshold</p> <p>TDC, QDC calc enable</p> <p>rise edge</p> <p>Set above noise</p> </div>

<p>timing type</p>	<p>Select the waveform for time stamp from CFD waveform and LED (raw waveform). Leading-Edge Timing (LET) It is the timing when a certain trigger level t is reached. Trigger acquisition timing is different if the wave height changes as a 'and b'.</p>
--------------------	--



Pic. 6: Concept of Leading-Edge Timing (LET)

Constant Fraction Discriminator Timing (CFD)











Pic. 7: Concept of Constant Fraction Discriminator Timing

Generate the waveforms "c, d", "e, f" and "g, h" in the table below for the different waveforms "a, b" in the figure above.


Waveforms c and d	Waveforms a and b multiplied by the CFD function and inverted
Waveforms e and f	Waveforms obtained by delaying the waveforms a and b by the CFD delay
Waveforms g and h	Waveforms obtained by adding waveforms c and e and waveforms obtained by adding d and f

CFD, which is the zero-cross timing of waveforms g and h, is characterized in that its value remains constant even if the wave height changes, provided that the rise time of the waveform is the same.

<p>CFD function</p>	<p>Magnification to reduce the original waveform for CFD waveform shaping.</p> <p>Set from 0.03 times, 0.06 times, 0.09 times, 0.12 times, 0.15 times, 0.18 times, 0.21 times, 0.25 times, 0.28 times, 0.31 times, 0.34 times, 0.37 times, 0.40 times, 0.43 times, 0.46 times.</p>  <p>The diagram shows a 'Raw Waveform' on the left. An arrow labeled 'x0.06' points to a smaller version of the waveform. A second arrow labeled 'x0.46' points to a larger version of the waveform. Blue double-headed arrows indicate the vertical scaling of the peaks.</p>
<p>CFD delay</p>	<p>Set the CFD delay time.</p> <p>This board sets in 1 ns step from 1 ns to 24 ns.</p>  <p>The diagram shows a 'Raw Waveform' on the left. An arrow labeled '2ns' points to a waveform with a small blue double-headed arrow indicating a narrow delay window. A second arrow labeled '24ns' points to a waveform with a larger blue double-headed arrow indicating a wider delay window.</p>
<p>CFD walk</p>	<p>Set the threshold for timestamping. The unit is digit.</p> <p>While viewing the waveform of "CFD" in wave mode, set the value near the zero-cross position.</p>  <p>The diagram shows a waveform with a 'Time Stamp' marker pointing to a peak. A horizontal blue line labeled 'WALK' is drawn across the waveform, with a small circle at its intersection with the zero-crossing point of the signal.</p>
<p>QDC sum or peak</p>	<p>Select the output format of QDC data.</p> <p>Select from PEAK value and SUM value.</p> <p>When PEAK is selected, the value of PEAK for the raw waveform is output as a QDC value</p> <p>When SUM is selected, FILTER is applied to the raw waveform and integrated value is output as QDC value</p>  <p>The diagram shows a 'Raw Waveform' on the left. An arrow labeled 'Filter' points to a smoothed, bell-shaped curve on the right. A blue shaded area under the curve represents the integrated value. Arrows point from the text above to the peak of the raw waveform and the shaded area of the filtered curve.</p>


<p>QDC pre trigger</p>	<p>Set the timing to start waveform shaping for integral value calculation. Select from 0 ns, -8 ns, -16 ns, -32 ns, -40 ns.</p> <p>Start integration from the previous time by the pretrigger setting value</p> 
<p>QDC filter</p>	<p>Set the time constant for shaping the waveform for integral value calculation. The setting is selected from Ext, 10 ns, 20 ns, 50 ns, 100 ns and 200 ns.</p> 
<p>QDC integral range</p>	<p>Select QDC integration time. The range is 0 ns to 32000 ns.</p> 
<p>QDC full scale</p>	<p>Set the gain of QDC data. The setting is selected from 1/1, 1/2, 1/4, 1/8, 1/16, 1/32, 1/64, 1/128, 1/256, 1/512. Make the QDC value be less than 8191.</p> 
<p>QDC LLD</p>	<p>Set the Lower Level Discriminator (LLD) of QDC. The unit is digit. Integration values below this threshold do not acquire timestamp data and integration value data. Set to a value smaller than ULD. The setting range is 0 to 8191.</p>
<p>QDC ULD</p>	<p>Set the Upper Level Discriminator (ULD) of QDC. The unit is digit. Integrated values above this threshold do not acquire time stamp data and integrated value data. Set to a value larger than LLD. The setting range is 0 to 8191.</p>
<p>OR eable</p>	<p>Set the OR output to the front panel AUX terminal. It can be set for each CH. One pulse TTL logic is output for the event after passing LLD and ULD. (optional) * Not shown in the photo on page 13.</p>


OR length	Set the TTL logic pulse width. It can be set from 8ns to 1000ns. (optional) * Not shown in the photo on page 13.
-----------	---



Pic. 5-2: config tab, lower left part

mode	Select from hist, list, wave, and list_com
	hist: Integrate the input signal and display the spectrum.
	wave: Digitize the input signal and display the waveform.
	list: Regarding the input signal, time information, CH information and integration information are regarded as one event, and output and stored as a binary file. Also used when acquiring the time spectrum.
	list_com: Use when measuring in time between boards. Use CH1 as a common signal input terminal and input a fast-rising pulse with little jitter after start.
measurement mode	Select from real time and live time. The measurement is completed in the selected time mode.
measurement time	Set the measurement time. It is up to 8760 hours.
list read byte	Set the unit reading number. It will be fixed at 16,000 bytes.
time spectrum on/off	Select whether to display time spectrum while acquiring list data in list mode. If you want to obtain only list data, uncheck it. NOTE: Selecting that check box during high count measurement will delay the acquisition of list data.
energy spectrum on/off	Select the presence or absence of spectrum display while acquiring list data in list mode. If you want to obtain only list data, do not select the check box. NOTE: Selecting that check box during high count measurement will delay the acquisition of list data.






Pic. 6: Setting example of rise start cnt and rise stop cnt

fall start cnt	<p>Start position of the target range of the integrated value FALL at the falling edge.</p> <p>Set the start position of integration range from the position exceeding threshold.</p> <p>The setting range is 1 to 16383 (16383 ns = 16383 × 1 ns).</p> <p>Set a value smaller than "fall stop cnt" described below</p>
fall stop cnt	<p>End position of the target range of integral value FALL of falling part.</p> <p>Set the range of integration from the above "fall start cnt".</p> <p>The setting range is 1 to 16383 (16383 ns = 16383 × 1 ns).</p> <p>Set a value larger than the "fall start cnt" described above.</p>

Example of calculation of FALL value:

In the case of set threshold: 50, fall start cnt: 5, fall stop cnt: 25, PSA full scale: In the case of 1/1, the FALL value exceeds threshold and integrates the blue frame part in the figure below for the 5th to 25 points. You The integral value is multiplied by PSA full scale to obtain the FALL value of the list data.




Pic. 7: Setting example of fall start cnt and fall stop cnt

total start cnt	Start position of the target range of the waveform integral value TOTAL. Set the range in front of the position beyond the threshold. The setting range is 1 to 498 (498 ns = 498 × 1 ns).
total stop cnt	End position of the target range of the waveform integral value TOTAL. Set the range for integration from the above "total start cnt". The setting range is 1 to 16383 (16383 ns = 16383 × 1 ns).

Example of calculation of TOTAL value:


In the case of set threshold: 50, total start cnt: 5, total stop cnt: 50, PSA full scale: In case of 1/1, integrate the red frame part in the figure below 50 points from 5 points before the position of threshold . The integral value is multiplied by PSA full scale to obtain the TOTAL value of the list data.



Pic. 8: Setting example of total start cnt and total stop cnt

PSA full scale	Set the reduction ratio of RISE value, FALL value and TOTAL value of list data
----------------	--

5. 3. file tab




Pic. 9: file tab

Setting section for saving.

histogram save	Save the histogram data displayed in the "spectrum tab" at the end of measurement to a file. The file save destination will be in the format described later. It is valid only when "hist" is selected in "mode"
histogram continuous save	Set whether to save histogram data to file continuously at set time interval. It is valid only when "hist" is selected in "mode".
histogram file path	Set absolute path of histogram data file. No extension is also possible. NOTE: It is not saved with this file name, but based on this file name, it will be the following format. Example: When "C: ¥ Data ¥ histogram.csv" is set in "histogram file path" and "10" is set in "histogram file save time (sec)", and the date and time is 2010/09/01 12:00:00 Will start saving data with the file name "C: ¥ Data ¥ histogram_20100901_120000.csv". After 10 seconds, save the file as "C: ¥ Data ¥ histogram_2010001_1200.csv". * The above "120010" may be "120009" or "120011".
hisotogram file save time (sec)	Set the time interval for continuous storage of histogram data. The unit is seconds. The setting range is 5 seconds to 3600 seconds.
list save	Set whether to save list data to file. Valid only when "list" is selected in "mode" in the Config tab.
list file number	Set the start number of the number added to the list data file. From 0 to 99,999. It will be reset to 0 if it exceeds 999999.

list file size (Byte)	Set the maximum file size of the list data file. If this size is exceeded while saving list data, the file is closed, and data saving is continued with a new file name that is incremented by one "list file number". "File size (Byte)" located on the right side of the setting displays the size of the file currently being saved.
list header ON/OFF	Set header ON / OFF when acquiring list data. The header is an IP address. When OFF, data without IP address header is saved.




5. 4. wave tab




Pic. 10: wave tab

Settings related to waveform display.

Graph	Waveform graph. When "wave" is selected in "mode" in the "config" tab, the waveform is displayed.
on/off	Set whether to display the waveform
CH	Select the CH of the waveform to be displayed.
Type	Select the type of waveform to be displayed. raw : Waveform digitized by the ADC and BLR processed CFD : Waveform shaped waveform of CFD. Filter : Waveform integrated by QDC PTG : Square wave at piled up timing
trigger edge	Select the trigger polarity. In general, select pos.
threshold	Set the trigger threshold. *You can also set the cursor in the graph.
trigger point	Sets the waveform display start point. *You can also set the cursor in the graph.
trigger SIG	Select a SIG (Signal) to trigger. Normally, select SIG1.
wave compress	Set the time scale compression of X axis. Used to display waveforms with long fall times.
wave free run	When the check box is not selected, the triggered waveform is displayed, and when it is checked, the trigger free waveform is displayed. It can also be used to view baseline levels and noise levels.


accumulation	Selection of valid / invalid of superposition of waveform data
XY Scale	Button to adjust the scale of X axis and Y axis. Expansion is + (plus), reduction is-(minus)
X axis calibration	Select the unit of X axis
Y axis calibration	Select the unit of Y axis. *Please use the mV display as a reference.
X-axis range	When you right-click on the X-axis and select "Auto scale", it becomes auto scale. When deselected, it will not be in automatic scale, and the minimum value and maximum value of X axis will be fixed. If you want to change the minimum or maximum value, you can place the mouse pointer over the number you want to change and click or double-click to change it.
Y-axis range	When you right-click on Y-axis and select "Auto scale", it becomes auto scale. When deselected, it will not be in automatic scale, and the minimum and maximum values of the Y axis will be fixed. If you want to change the minimum or maximum value, you can place the mouse pointer over the number you want to change and click or double-click to change it.
	Cursor movement tool. The cursor can be moved on the graph when setting the ROI.
	Zoom. When selected, the following six zooms in and zoom out options can be selected and executed.
	Pan tool. You can grab the plot and move it over the graph.



Pic. 11: For zooming in and out

1	Rectangle	Use this option to click the point on the display that you want to be the corner of the zoom area and drag the tool until the rectangle occupies the zoom area.
2	X-Zoom	Zooms in to the area of the graph along the X axis.
3	Y-Zoom	Zooms in to the area of the graph along the Y axis.
4	Fit Zoom	Automatically scales all X and Y scales on the graph.
5	Zoom out around point	Click the center point to zoom out.
6	Zoom in around point	Click the center point to zoom in.

5. 5. spectrum tab




Pic. 12: spectrum tab

Settings related to spectrum display.

Graph	Energy spectrum. The spectrum is displayed when "hist" is selected for "mode" in the "config" tab or when "list" is selected for "mode" in the "config" tab and "spectrum ON / OFF" is enabled.
Check box	Sets whether to display the histogram for each channel on the graph.
ROI CH	Select the CH number to apply ROI (Region Of Interest). Up to 8 (eight) ROIs can be set for one CH signal
ROI start (ch)	Set the start position of ROI. The unit is ch
ROI end (ch)	Set the end position of the ROI. The unit is ch
energy	Defines the energy value of peak position (ch). In the case of 60 Co, it is set to 1173 (keV) or 1332 (keV). When "ch" is selected in "calibration", the peak between the ROIs is detected, keV / ch is calculated from the peak position (ch) and the set energy value, and the calculated half width is calculated.
calibration	Select the unit of X axis. The label of the X axis is also changed along with the setting. ch : Unit ch (channel) display. Units such as "FWTM" of "FWTM" of ROI are arbitrary. eV : Unit eV display. Calculate the slope a and intercept b of the linear function $y = ax + b$ so that ch becomes eV by two-point calibration of two types of peaks (center value) and energy value in one histogram, and set as the X axis You The unit such as "FWTM" of "FWTM" of ROI is "eV".

	<p>keV: Unit keV display. Calculate the slope a and intercept b of the linear function $y = ax + b$ so that ch becomes keV by two-point calibration of two kinds of peaks (center value) and energy value in one histogram, and set as the X axis You The unit such as "FWTM" of "FWTM" of ROI is "keV". Example: When there are 1173.24keV of 60Co in 5717.9ch and 1332.5keV of 60Co in 6498.7ch, a is calculated automatically as 0.20397 and b is 6.958297 from 2-point calibration.</p> <p>manual: Set the slope a and intercept b of the linear function $y = ax + b$, and the unit label to the X axis. The unit is set arbitrarily.</p>
Y mapping	<p>Select mapping of Y axis of graph. The label of the Y axis is also changed along with the setting. linear: straight line log: log</p>
smoothing	<p>Smoothing (5-point second order smoothing) function to calculate the half-width when there are few statistics.</p>
simple count view	<p>You can easily read the count displayed on the graph.</p>
gauss fit	<p>Apply Gaussian fitting to the spectrum.</p>

5. 6. timespectrum tab



Pic. 2: timespectrum tab

Settings related to timespectrum display. It is a measurement limitation in the board.

*Timespectrum is generated based on the list data acquired in list mode.


Graph	Time difference spectrum. If "list" is selected for "mode" in the "config" tab and "timespectrum on / off" is selected, the time difference spectrum is displayed during measurement.
Check box	Select whether to display spectrum or not.
Config part	Setting of time spectrum. start CH: Select the CH number for acquiring the start timing. stop CH: Select the CH number for acquiring stop timing. gain: You can select from 1x to 1 / 128x. At 1-time, full scale about 780 ns (about 3.9 ps per 1 digit), 1/128 full scale is about 100 μs (0.5 ns per 1 digit). coinc offset: Sets the offset in units of 1 ns. coinc time: Set the coincidence time in units of 1 ns. If the time difference between event detection in the above "start CH" and "stop CH" is within this setting range, it is considered as coincidence (simultaneous) and it is considered as valid data.
ROI	Settings related to calculation. ROI START: ROI start channel ROI END: end channel of ROI FWHM: The calculated half width is displayed. FWTM: Displays the calculated full value range.
Xscale	Select the unit of X axis, "ch" channel or "ns" display.

6. Measurement

As an example, the operation procedure of energy spectrum measurement, list measurement, and time spectrum measurement when using LaBr3 (Ce) detector (hereinafter referred to as detector) is described.

6. 1. Energy Spectrum measurement

(1) Environment




Pic. 3: Configuration for energy spectrum measurement

	Check that all devices (VME rack, High voltage power supply and PC) are off.
	Connect the detector and HV with the SHV connector cable.
	Connect the anode output signal from the detector to CH1 of APV8108-14 with LEMO connector coaxial cable. In the case of a BNC connector, use a BNC-LEMO conversion adapter.
	Connect APV8108-14 and PC with LAN cable.
	Power on the VME rack.
	Turn on the PC. Launch the application.
	Turn on the high voltage power supply and apply the voltage according to the detector.
	In this example, Cs-137 is used for the radiation source.


(2) Waveform measurement

First, check the signal from the detector that is input in waveform mode. Make the following settings in the "config" tab, and then click the "Config" menu.



Pic.4 Setting of waveform measurement

Open the "wave" tab, check the settings shown below, click "Clear" in the menu and then click "Start". The waveform from the detector is displayed on the graph.




Pic.5: Screen of waveform measurement

Note the following points

If no signal is displayed, the trigger may not have been applied.

Select "wave free run" in the "wave" tab to check the baseline, click "Config" on the menu bar, then "Clear" and finally "Start".

You can see how many wave heights are at the baseline and roughly.



Pic. 6: Screen of baseline check

Next, deselect "wave free run", gradually increase the "threshold" value from about 10, and confirm the "threshold" value at which the waveform can be captured firmly as shown on the previous page. This value will be used for later configuration.


* Check if the wave height is too large to saturate.

Decrease the applied high voltage to reduce the amplitude of the input signal to the device.

Measured data can be saved with "save wave" in "File" of the menu bar.


(3) Measurement of Energy Spectrum

For spectrum measurement, click the "Config" on the menu bar after making the following settings on the "config" tab. Set the "threshold" value set during waveform measurement to "threshold" in the "config" tab.



Pic. 7: Config tab

Open the "spectrum" tab, check the settings shown below, then click "Clear" and then "Start" in the menu bar. The following spectrum is displayed after execution.



Pic. 8: Screen of Energy spectrum measurement

Please note the following points.

- Select "spectrum on / off" CH1 and display the spectrum of CH1.
- When analyzing peaks, set the ROI. Refer to "5.5. Spectrum tab" for details.

Measured data is saved with "save histogram" in "File" of the menu bar.

When you finish measurement, click "Stop" on the menu bar.

6. 2. List measurement

(1) Environment




Fig. 9: Configuration for energy spectrum measurement

	Check that all devices (VME rack, High voltage power supply and PC) are off.
	Connect the detector and HV with the SHV connector cable.
	Connect the anode output signal from the detector to CH1 of APV8108-14 with LEMO connector coaxial cable. In the case of a BNC connector, use a BNC-LEMO conversion adapter.
	Connect APV8108-14 and PC with LAN cable.
	Power on the VME rack.
	Turn on the PC. Launch the application.
	Turn on the high voltage power supply and apply the voltage according to the detector.
	In this example, Cs-137 is used for the radiation source.

(2) Confirmation of input waveform

Confirm the same as in “6. 1. Energy Spectrum Measurement (2) Waveform Measurement” above.

(3) Confirmation of Energy spectrum

Confirm the same as in “6. 1. Energy Spectrum Measurement (3) Measurement of Energy Spectrum” above.

NOTE: the following points in this software.

output rate (cps)	It is the number of events to be acquired in one second, and it is checked whether it is too low or too high for the assumption (See No. 1 in the figure on the next page). In list mode, 16 bytes of data are acquired for each event. For example, if "output rate (cps)" is 500 kcps, data of 8 MB / s (500 kcps x 16 Bytes) will be saved per second.
spectrum tab	Check that there is no abnormality in the shape of the spectrum, and if excessive noise data is acquired. (See No. 2 on the next page)




Fig. 10: Point to make note of before measurement in list mode

(4) List measurement

Start list measurement. Set "mode" to "list" in the "config" tab.




Fig. 11: config tab

To save list data, set the following items in the "file" tab.

list save	Check
list file path	File path to be a reference
list file number	Set arbitrarily in the range from 0 to 999999. Please be careful not to duplicate.
list file size (Byte)	Size of list data file. When this size is exceeded, "list file number" is automatically incremented by one and saved to a new file.




Fig.12: Settings for saving list data in file tab

Click in the order of "Config" "Clear" "Start" on the menu bar.

After execution, when the event is detected and list data is acquired, the following "file size (Byte)" is increased.




Fig. 13: Screen of list data measurement and saving

When you finish measurement, click "Stop" on the menu bar.

6. 3. Time Spectrum measurement

(1) Environment




Fig. 14: Configuration of time spectrum measurement

	Check that all devices (VME rack, High voltage power supply and PC) are off.
	Connect the detector and HV with the SHV connector cable.
	Connect the anode output signal from the detector to CH1 of APV8108-14 with LEMO connector coaxial cable. In the case of a BNC connector, use a BNC-LEMO conversion adapter.
	Connect APV8108-14 and PC with LAN cable.
	Power on the VME rack.
	Turn on the PC. Launch the application.
	Turn on the high voltage power supply and apply the voltage according to the detector.
	In this example, Na-22 is used for the radiation source.

(2) Waveform measurement

Confirm the same as in “6. 1. Energy Spectrum Measurement (2) Waveform Measurement” above.

(3) Energy Spectrum measurement

While checking the signal from the detector, specify the range of energy for time measurement.

First, perform the energy spectrum measurement with the following settings. Make the following settings on the "config" tab, and then click "Config" on the menu bar.




Fig. 15: Energy spectrum measurement setting before time spectrum measurement (full energy range)

Open the "spectrum" tab and click "Clear" and then "Start" in the menu bar.

The following spectrum is displayed after execution. While checking the shape of the spectrum and counting, set the target of the peak range with "ROI start" and "ROI end".




Fig. 16: Energy spectrum measurement setting before time spectrum measurement (full energy range)

Next, set the following to narrow down to the energy to be measured (in this case, the 511 keV peak of Na-22).

Set the reference values for "ROI start" and "ROI end" on the previous page in the "config" tab in the red frame below. Enter a value for "ROI start" for "QDC LLD" and a value for "ROI end" for "QDC ULD."




Fig. 17: Energy spectrum measurement setting before time spectrum measurement (Setting of energy range narrowing)

Open the "spectrum" tab and click "Clear" and then "Start" in the menu bar. The following spectrum is displayed after execution. The energy peaks shown in the figure below, which are narrowed down in the range of "QDC LLD" and "QDC ULD", are displayed.




Fig. 18: Energy spectrum measurement setting before time spectrum measurement (Setting of energy range narrowing)

(4) Time Spectrum measurement

If you want to measure the spectrum, select "timespectrum ON / OFF", make the following settings on the "config" tab, and then click "Config" on the menu bar.

NOTE: When "mode" is selected in "list" mode. Note that if you use high counting in this mode, the load on your computer may be unstable and the behavior may become unstable.




Fig. 19: Setting of time spectrum measurement




Fig. 20: Setting of time spectrum measurement

Open the "spectrum" tab and click "Clear" and then "Start" in the menu bar. The following spectrum is displayed after execution. The time resolution "FWHM (ps)" is calculated by setting the "ROI" section at the lower right of the screen.




Fig.21: Time spectrum measurement




Fig. 22: Time spectrum measurement (When the horizontal axis is magnified)

When you finish measurement, click "Stop" on the menu bar.

7. File

7. 1. Histogram data file

(1) File format

CSV text format, separated by commas

(2) File name

Set arbitrarily

(3) Component

It consists of 4 parts: "Header" part, "Calculation" part, "Status" part and "Data" part

Header part

Measurement mode	Measurement mode.
Measurement time	Measurement set time. Unit is seconds
Real time	Real time
Start Time	Measurement start time
End Time	Measurement end time
* Save for each CH	
POL	Polarity
TGE	Waveform display trigger CH
TGC	Waveform acquisition polarity
RJT	Waveform acquisition threshold
CCF	CFD Function
CDL	CFD Delay
CWK	CFD walk
CTH	CFD threshold
FLK	Baseline time constant
PTS	QDC pre-trigger
LIG	QDC filter time constant
LIT	QDC thumb or peak
AFS	QDC integral reduction
CLD	QDC Lower Limit Discriminator
CUD	QDC Upper Limit Discriminator
TTY	Timing type
MOD	mode
MTM	Measurement time
MEMO	Memo or Notes

Calculation part

*The following is saved for each ROI

ROI_ch	Input channel number targeted for ROI.
ROI_start	ROI start position (ch)
ROI_end	ROI end position (ch)
energy (keV)	ROI setting energy (keV)
peak (ch)	Peak position between ROIs (ch)
centroid (ch)	Center position between ROIs (ch)
peak (count)	Peak ch count between ROIs
gross (count)	Total number of counts between ROIs
gross (cps)	Cps of counts between ROIs
net (count)	Total number of counts minus background between ROIs
net (cps)	Sum of cps of the number of counts minus background between ROIs
FWHM (ch)	Half width between ROIs (ch)
FWHM (%)	Resolution between ROIs (%)
FWHM (keV)	Half-width between ROIs (keV)
FWTM (keV)	Full width between ROIs (keV)

Status part

*The following is saved for each ROI

output count	Output count
output rate	Output count rate
dead time	Dead time ratio

Data part

Histogram data for each channel. Up to 8192 points.

7. 2. Wave data file

(1) File format

CSV text format, separated by commas

(2) File name

Set arbitrarily

(3) Component

It consists of 4 parts: "Header" part, "Calculation" part, "Status" part and "Data" part

Header part

Measurement mode	Measurement mode.
Measurement time	Measurement set time. Unit is seconds
Real time	Real time
Start Time	Measurement start time
End Time	Measurement end time
* Save for each CH	
POL	Polarity
TGE	Waveform display trigger CH
TGC	Waveform acquisition polarity
RJT	Waveform acquisition threshold
CCF	CFD Function
CDL	CFD Delay
CWK	CFD walk
CTH	CFD threshold
FLK	Baseline time constant
PTS	QDC pre-trigger
LIG	QDC filter time constant
LIT	QDC thumb or peak
AFS	QDC integral reduction
CLD	QDC Lower Limit Discriminator
CUD	QDC Upper Limit Discriminator
TTY	Timing type
MOD	mode
MTM	Measurement time
MEMO	Memo or Notes

Status part

*The following is saved for each ROI

output count	Output count
output rate	Output count rate
dead time	Dead time ratio

Data part

Waveform data of device being displayed

7. 3. List data file

(1) File format

Binary, network byte order (big endian, MSB first)

(2) Configuration

The APV8108-14 sends binary data of the following format to the PC sequentially in the list mode.

Bit127		112	
TOTAL[15..0]			
111		96	
FALL[15..0]			
95		80	
RISE[15..0]			
79		64	
TDC[55..40]			
63		48	
TDC[39..24]			
47		32	
TDC[23..8]			
31		24	16
TDC[7..0]		TDCFP[7..0]	
15	13	12	0
CH[2..0]		QDC [12..0]	

Fig. 23: list data format (16 Byte (128-bit))

Details of list data

Bit 127 to Bit 112	TOTAL (total waveform integral) value. Unsigned 16-bit integer.
Bit 111 to Bit 96	FALL (falling waveform partial integration) value. Unsigned 16-bit integer.
Bit 95 to Bit 80	RISE (Rise of waveform partial integration) value. Unsigned 16-bit integer.
Bit 79 to Bit 24	TDC count. 56-bit. 1 ns per bit.
Bit 23 to Bit 16	TDCFP (fractional part) count. 8-bit. 3.90625 ps per bit. Interpolation between sampling points (1 ns 256 256 = 3.90625 ps)
Bit 15 to Bit 13	CH number. 0: CH1, 1: CH2, 2: CH3, 3: CH4, 4: CH5, 5: CH6, 6: CH7, 7: CH8
Bit 12 to Bit 0	QDC integrated value. Unsigned 13-bit integer. The collected waveform is filtered, and the integrated value of the waveform between the set range from where the threshold is exceeded.

8. Command

8. 1 Overview

Configuration and data acquisition for the APV 8108 are performed via TCP / IP and UDP via Ethernet. Since special libraries are not used, DPP can be controlled by any application if it conforms to the communication format (command).


A communication board is mounted on DPP.

Model Number	Standard	Protocol	Command method
APG5107	1000 Mbps	TCP/IP and UDP	Address + parameter

This chapter describes commands when the APG5107 is mounted as a communication board.

The APG5107 uses SiTCP to realize high-speed data communication. SiTCP is a technology for connecting devices developed by the High Energy Accelerator Research Organization (<http://www.kek.jp/ja/>) (hereinafter referred to as KEK), which is a university shared use corporation, to Ethernet. The technology has been transferred to Bee Beans Technologies Co., Ltd. (<http://www.bbtech.co.jp>, hereinafter BBT), which is a venture company from KEK. We use SiTCP under license from BBT. For details on SiTCP and data transmission / reception, please refer to each manual on the BBT company website.

The command types are roughly divided into two types: "Configuration and Status" and "Data". In SiTCP, two protocols, TCP / IP and UDP, are operating so that these two types of commands can be sent and received without competition, and each has a communication port on the device side defined. Configuration and Status are UDP and port number is 4660 by default. Data is TCP / IP and the port number is 24 by default.



The format and type of command are described below.

8. 2 Format of command

The command format is divided into the case of Configuration write, the case of Status read, and the case of Data read. Each consists of "header", "address", "parameter" and "data".

"Header part" contains 6 items of Ver / Type / CMD / FLAG / ID and Data Length conforming to the SiTCP specification. In DPP, Data Length (data length) is fixed 2-byte and the size of the header part is 4-byte.

"Address part" is the 4-byte address of the DPP register.

"Parameter part" is a 2-byte value set in the DPP register.

"Data part" is measurement data from DPP.

8. 3 Type of command

(1) Config command




Fig. 1: In case of Config command

The Config command responds from the DPP to the configuration from the PC.

【Setting from PC】

"Header part" is 4-byte, hexadecimal number "FF800702".

The contents are F (Ver.) F (Type) 8 (CMD) 0 (FLG) 07 (ID) 02 (Data Length).

"Address part" sets 4-byte, parameter address value.

"Parameter part" is 2-byte, the parameter value to be set.

【Response from DPP】

"Header part" is 4-byte, hexadecimal number "FF880702".

The contents are F (Ver.) F (Type) 8 (CMD) 8 (FLG) 07 (ID) 02 (Data Length).

If normal, the ACK bit of FLG becomes 1 and becomes 8.

"Address part" returns 4 bytes, the address value of the set parameter.

(2) **Status command**




Fig.2: In case of Status command

The Status command responds to requests from the PC from DPP.

【Requirements from PC】

"Header part" is 4-byte, hexadecimal number "FFC00602".

The contents are F (Ver.) F (Type) C (CMD) 0 (FLG) 06 (ID) 02 (Data Length).

"Address part" sets 4-byte, the address value of status data.

【Response from DPP】

"Header part" is 4-byte, hexadecimal number "FFC 80602".

The contents are F (Ver.) F (Type) C (CMD) 8 (FLG) 06 (ID) 02 (Data Length).

If normal, the ACK bit of FLG becomes 1 and becomes 8.

"Address section", 4-byte, the address value of the requested parameter is returned.

"Data part" is 2-byte, status data value.

(3) **Data command**

Read list data to PC with Data command.

First, list data is requested by UDP, and DPP responds by TCP / IP.

【TCP / IP DPP data transmission】

When the list mode is set and started by UDP, the list data starts to accumulate on the communication board.

【TCP / IP PC data reception】

It is possible to read only arbitrary bytes by TCP.

Since list data is in units of 16 bytes, read in units of 16 bytes.

8. 4 List of command

No.	Type	No. of port	Address (Hexadecimal)	Content	Setting range (Digit)	Operation	Command length (Byte)	Response (Byte)
1	CH setting	UDP 4660	B40001DE	Input waveform type	0, 1	Configuration	10	10
						Setting request	8	10
2			B400011A	Input polarity switching	0, 1	Configuration	10	10
						Setting request	8	10
3			B4000160	CFD Function	1..15	Configuration	10	10
						Setting request	8	10
4			B4000162	CFD Delay	0..23	Configuration	10	10
						Setting request	8	10
5			B4000164	CFD Walk	0.. 2 ¹⁰ -1	Configuration	10	10
						Setting request	8	10
6			B4000166	Threshold	0.. 2 ¹³ -1	Configuration	10	10
						Setting request	8	10
7			B400016E	Baseline restoration filter	0..254	Configuration	10	10
						Setting request	8	10
8			B40001C0	QDC pre-trigger	0..4	Configuration	10	10
	Setting request	8				10		
9	B40001C6	QDC Filter	0..5	Configuration	10	10		
				Setting request	8	10		
10	B40001C8	QDC sum/peak	0, 1	Configuration	10	10		
				Setting request	8	10		
11	B400010C	QDC full scale	0..9	Configuration	10	10		
				Setting request	8	10		
12	B40001DC	QDC integral range	1.. 2 ¹² -1	Configuration	10	10		
				Setting request	8	10		
13	B4000168	QDC LLD	0.. 2 ¹³ -1	Configuration	10	10		
				Setting request	8	10		
14	B400016A	QDC ULD	0.. 2 ¹³ -1	Configuration	10	10		
				Setting request	8	10		
15	B40001D0	Time stamp timing	0, 1	Configuration	10	10		
				Setting request	8	10		

No.	Type	No. of port	Address (Hexadecimal)	Content	Setting range (Digit)	Operation	Command length (Byte)	Response (Byte)
16	CH setting	UDP 4660	B40001D8	PSA falling start position	1..16383	Configuration	10	10
						Setting request	8	10
17			B40001DA	PSA falling end position	1..16383	Configuration	10	10
						Setting request	8	10
18			B40001E8	PSA rising start position	1..498	Configuration	10	10
						Setting request	8	10
19			B40001EA	PSA rising end position	1..16383	Configuration	10	10
						Setting request	8	10
20			B40001EC	PSA overall start position	1..498	Configuration	10	10
						Setting request	8	10
21	B40001EE	PSA overall end position	1..16383	Configuration	10	10		
				Setting request	8	10		
22	B40001D6	PSA reduction ratio	0..9	Configuration	10	10		
				Setting request	8	10		
23	B4000176	Input delay	0..511	Configuration	10	10		
				Setting request	8	10		

NOTE: The above address is that of CH1. The start address of CH1 / CH5 is B4000100 / B4008100. The start address of CH2 / CH6 is B4000200 / B4008200. The address to which 0x100 is added in this way becomes the top of each CH setting.

No.	Type	No. of port	Address (Hexadecimal)	Content	Setting range (Digit)	Operation	Command length (Byte)	Response (Byte)				
24	Single setting	UDP 4660	B4004000	Mode	0,1,2,5	Configuration	10	10				
						Setting request	8	10				
25			B4004002	Measurement mode	0,1	Configuration	10	10				
						Setting request	8	10				
26			B4004006 B4004008 B400400A B400400C	Measurement time setting	0..2 ⁵⁴ -1	Configuration	10*4 times	10*4 times				
						Setting request	8*4 times	10*4 times				
27			B4004004	Start measurement	0, 1	Configuration	10	10				
28			B4004090	Time and data clear	0, 1	Configuration	10	10				
29			B4004028	Time clear	0, 1	Configuration	10	10				
30	Status	UDP 4660	B4000004	Measurement status	-	Status request	10	10				
31 *			B4000120 B4000122	Output count total			10*2 times	10*2 times				
							32 *	B4000130 B4000132	Output count rate	10*2 times	10*2 times	
33 *			B4000144 B4000146 B4000148 B400014A	Live count						10*3 times	10*3 times	
							34 *	B40001E0 B40001E2 B40001E4 B40001E6	Dead count	10*4 times	10*4 times	
										35	B400000E B4000010 B4000012 B4000014	Real time
36			Configuration	UDP 4660			B400009A B400809A	Histogram CH date	0..7	Configuration	10	32768

NOTE: The address marked with * in the above number column is for CH1. The addresses for CH1 to CH4 address 0x100 and CH6 to 8 CH5 start address: 4008100 plus 0x100 are the values for each CH setting. (Legend CH1: B4000146 CH CH2: B4000246, ..., CH4: B4000446, CH5: B4008146, .. CH8: B4008446)

8. 5 Explanation of command

Setting of CH

(1) Select input waveform type

Description	Select input waveform type
Address	0xB40001DE
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 or 1 0: Normal waveform * Factory default 1: NIM waveform

(2) Select Input polarity switching

Description	Select input polarity switching
Address	0xB40001A
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 or 1 0: Negative polarity * Factory default 1: Positive polarity

(3) CFD Function

Description	CFD function settings. Signal reduction factor used for CFD waveform calculation.
Address	0xB4000160
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 15 1: 0.03 times, 2: 0.06 times, 3: 0.09 times, 4: 0.12 times, 5: 0.15 times, 6: 0.18 times, 7: 0.21 times*, 8: 0.25 times, 9: 0.28 times, 10: 0.31 times, 11: 0.34 times, 12: 0.37 times, 13: 0.40 times, 14: 0.43 times, 15: 0.46 times

(4) CFD Delay

Description	CFD delay settings. Delay time of inverted signal used for CFD waveform calculation.
Address	0xB4000162
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 23 0: 1 ns, 1: 2 ns, 2: 3 ns, 3: 4 ns, 4: 5 ns*, 5: 6 ns, 6: 7 ns, 7: 8 ns, 8: 9 ns, 9: 10 ns, 10: 11 ns, 11: 12 ns, 12: 13 ns, 13: 14 ns, 14: 15 ns, 15: 16 ns, 16: 17 ns, 17: 18 ns, 18: 19 ns, 19: 20 ns, 20: 21 ns, 21: 22 ns, 22: 23 ns, 23: 24 ns

(5) CFD Walk

Description	Setting value for timestamping ADC data
Address	0xB4000164
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 1023 10 * Factory default

(6) Theshold

Description	Input waveform threshold setting
Address	0xB4000166
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 8191 100 * Factory default

(7) Baseline restorer filter

Description	Baseline restorer filter time constant
Address	0xB400016E
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 254 0: ext 64: fast 128: 4 μ s 250: 85 μ s 252: 129 μ s 254: 260 μ s*

(8) QDC pre-trigger

Description	Set the timing to start the addition
Address	0xB40001C0
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 4 0: 0 ns 1: Addition starts minus 8 ns before the timing of the threshold 2: Addition starts minus 16 ns before the timing of the threshold* 3: Addition starts minus 24 ns before the timing of the threshold 4: Addition starts minus 32 ns before the timing of the threshold

(9) QDC Filter

Description	Filter time constant of the original waveform to be integrated by QDC
Address	0xB40001C6
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 5 0: ext 1: 10 ns* 2: 20 ns 3: 50 ns 4: 100 ns 5: 200 ns

(10) QDC sum / peak

Description	QDC output data type
Address	0xB40001C8
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 or 1 0: Output the peak value of the peak waveform as QDC data 1: Output the integrated value of the waveform applied with sum Filter as QDC data *

(11) QDC full scale

Description	QDC data gain
Address	0xB400010C
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 9 0: 1/1 1: 1/2 times 2: 1/4 times* 3: 1/8 times 4: 1/16 times 5: 1/32 times 6: 1/64 times 7: 1/128 times

(1 2) QDC Integral range

Description	Integration time setting including QDC pre-trigger setting time
Address	0xB40001DC
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 4095 1: 1 digit corresponds to 8 ns. * 25 (200 ns) Factory default To 4095: 32760 ns

(1 3) QDC LLD

Description	LLD of integral value of QDC
Address	0xB4000168
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 8191 10 * Factory default

(1 4) QDC ULD

Description	ULD of integrated value of QDC
Address	0xB400016A
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 8191 8000 * Factory default

(1 5) Time stamp timing

Description	Select which waveform to use for the time stamp of time information.
Address	0xB40001D0
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 or 1 0: The timing set by CFDWALK for the CFD waveform

(1 6) PSA falling start position

Description	Set start position of fall integration target range in PSA operation
Address	0xB40001D8
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 16383 * 1 digit: 1 ns

(17) PSA falling end position

Description	Set the end position of fall integration target range in PSA operation
Address	0xB40001DA
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 16383 * 1 digit: 1 ns

(18) PSA rising start position

Description	Set the start position of rising integration target range in PSA operation
Address	0xB40001E8
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 498 * 1 digit: 1 ns

(19) PSA rising end position

Description	Set the end position of rising integration target range in PSA operation
Address	0xB40001EA
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 16383 * 1 digit: 1 ns

(20) Overall PSA start position

Description	Set the start position of the whole integration target range in PSA operation
Address	0xB40001EC
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 498 * 1 digit: 1 ns

(21) Overall PSA end position

Description	Set the end position of the whole integration target range in PSA operation
Address	0xB40001EC
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 16383 * 1 digit: 1 ns

(22) PSA reduction ratio

Description	Set the reduction ratio of each integration result of rising, falling, and whole
Address	0xB40001D6
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0 to 9 0: 1/1, 1: 1/2, 2: 1/4, 3: 1/8, 4: 1/16, 5: 1/32, 6: 1/64, 7: 1/128, 8: 1/256, 9: 1/512

(23) Input delay

Description	Set delay of input waveform
Address	0xB4000176
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	1 to 511 * 8 ns delay per digit

Single setting

(24) Mode

Description	Operation mode. Select hist (histogram) mode, wave (waveform) mode, or list (list) mode
Address	0xB4004000
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0, 1 or 2 0: Histo mode 1: Waveform mode * 2: List mode 5: List common mode In list common mode, the GATE / VETO signal is invalid for CH1 only.

(25) Measurement mode

Description	Select real time or live time
Address	0xB4004002
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Channel setting
Range	0, 1 or 2 0: Real time * 1: Live time

(26) Measurement time setting

Description	Measurement time
Address	0xB4004006 (MSB) 0xB4004008 0xB400400A 0xB400400C (LSB)
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Single setting
Range	0 to $2^{54} - 1$ $2^{54} - 1$ * Factory default 8 ns per bit. The maximum setting range is 40031 hours from $(2^{54}-1) \times 8$ ns

(27) Measurement starts

Description	Set start or stop of measurement
Address	0xB4004004
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Single setting
Range	0 or 1 0: Measurement stop * 1: Measurement start

(28) Time and data clear

Description	Clear setting of time and data
Address	0xB4004090
Setting	Command length 10 Byte, Response 10 Byte
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Single setting
Range	When clearing, set data sequentially in the order of 0, 1, 0

(29) Time clear

Description	Time clear
Address	0xB4004028
Setting	Command length 10 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Single setting
Range	When clearing, set data sequentially in the order of 0, 1, 0

Status

(30) Measurement status

Description	Confirm measurement status
Address	0xB4000004
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Status
Range	0 or 1 0: Measurement stopped 1: During measurement

(31) Output count total

Description	Total count number processed by signal within QDC LLD and QDC ULD
Address	0xB4000120 (MSB) 0xB4000122 (LSB)
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Status (CH)
Contents	Total counts processed by signal processing

(32) Output count rate

Description	1 second output count rate
Address	0xB4000130 (MSB) 0xB4000132 (LSB)
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Status (CH)
Contents	Signaled counts per second

(33) Live count

Description	Total live count (8 ns / count)
Address	0xB4000144 (MSB) 0xB4000146 0xB4000148 0xB400014A (LSB)
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Status (CH)
Contents	The time obtained by subtracting the dead time from the real time per CH (Conversion to time is converted by count value x 8ns)

(3.4) Dead count

Description	Total dead count (8 ns / count)
Address	0xB40001E0 (MSB) 0xB40001E2 0xB40001E4 0xB40001E6 (LSB)
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Status (CH)
Contents	Signal processing time operating within QDC LLD, QDC ULD (Conversion to time is converted by count value x 8ns)

(3.5) Real time

Description	Real time (8 ns / count)
Address	0xB400000E (MSB) 0xB4000010 0xB4000012 0xB4000014 (LSB)
Request	Command length 8 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Status (CH)
Contents	Elapsed time after about 2.1 seconds from the start of measurement. Signals of CH1 to 8 will not be received for about 2.1 seconds after the start of measurement. (Conversion to time is converted by count value x 8ns)

(3.6) Histogram CH data

Description	Request Histogram Data
Address	0xB400009A 0xB400809A
Setting	Command length 10 Byte, Response 10 Byte
Port number	4660 (UDP)
Type	Status (CH)
Contents	0xB400009A setting 0: CH1 histogram request 1: CH2 histogram request 2: CH3 histogram request 3: CH4 histogram request 0xB400809A setting 0: CH5 histogram request 1: CH6 histogram request 2: CH7 histogram request 3: CH8 histogram request NOTE: When you want to acquire the histogram data of CH1, if you set data 0x0000 to address 0xB400009A, there is a data response of 32768 bytes of histogram data of CH1 in TCP.

8. 6 Setting command at startup and config

The application for the APV 8108-14 is set from the PC to DPP by combining the setting command described in "8.4 Command List" and the APV 8108-14 unique setting command when starting up after power on and when operating the Config menu before starting measurement. to hold.

The following shows an example of setting command list at startup and at configuration (does not include DPP response). Please refer to it when creating your own application.

Except for the setting commands in "8.4 Commands" and "8.5 Commands", these commands are APV8108-14 specific setting commands, so please be careful not to delete or change them.

1) Startup setup command list example

UDP send	0xFF800702B4004000001	UDP send	0xFF800702B40081620009
UDP send	0xFF800702B40040060000	UDP send	0xFF800702B40082620009
UDP send	0xFF800702B40040080000	UDP send	0xFF800702B40083620009
UDP send	0xFF800702B400400A2540	UDP send	0xFF800702B40084620009
UDP send	0xFF800702B400400CBE40	UDP send	0xFF800702B40001640019
UDP send	0xFF800702B400402E0100	UDP send	0xFF800702B40002640019
UDP send	0xFF800702B40040300001	UDP send	0xFF800702B40003640019
UDP send	0xFF800702B400403200FF	UDP send	0xFF800702B40004640019
UDP send	0xFF800702B400008C0007	UDP send	0xFF800702B40081640019
UDP send	0xFF800702B400011A0001	UDP send	0xFF800702B40082640019
UDP send	0xFF800702B400021A0001	UDP send	0xFF800702B40083640019
UDP send	0xFF800702B400031A0001	UDP send	0xFF800702B40084640019
UDP send	0xFF800702B400041A0001	UDP send	0xFF800702B4000166001E
UDP send	0xFF800702B400811A0001	UDP send	0xFF800702B4000266001E
UDP send	0xFF800702B400821A0001	UDP send	0xFF800702B4000366001E
UDP send	0xFF800702B400831A0001	UDP send	0xFF800702B4000466001E
UDP send	0xFF800702B400841A0001	UDP send	0xFF800702B4008166001E
UDP send	0xFF800702B40040360000	UDP send	0xFF800702B4008266001E
UDP send	0xFF800702B400010C0004	UDP send	0xFF800702B4008366001E
UDP send	0xFF800702B400020C0004	UDP send	0xFF800702B4008466001E
UDP send	0xFF800702B400030C0004	UDP send	0xFF800702B4000168001E
UDP send	0xFF800702B400040C0004	UDP send	0xFF800702B4000268001E
UDP send	0xFF800702B400810C0004	UDP send	0xFF800702B4000368001E
UDP send	0xFF800702B400820C0004	UDP send	0xFF800702B4000468001E
UDP send	0xFF800702B400830C0004	UDP send	0xFF800702B4008168001E
UDP send	0xFF800702B400840C0004	UDP send	0xFF800702B4008268001E
UDP send	0xFF800702B40001600007	UDP send	0xFF800702B4008368001E
UDP send	0xFF800702B40002600007	UDP send	0xFF800702B4008468001E
UDP send	0xFF800702B40003600007	UDP send	0xFF800702B400016A1F40
UDP send	0xFF800702B40004600007	UDP send	0xFF800702B400026A1F40
UDP send	0xFF800702B40081600007	UDP send	0xFF800702B400036A1F40
UDP send	0xFF800702B40082600007	UDP send	0xFF800702B400046A1F40
UDP send	0xFF800702B40083600007	UDP send	0xFF800702B400816A1F40
UDP send	0xFF800702B40084600007	UDP send	0xFF800702B400826A1F40
UDP send	0xFF800702B40001620009	UDP send	0xFF800702B400836A1F40
UDP send	0xFF800702B40002620009	UDP send	0xFF800702B400846A1F40
UDP send	0xFF800702B40003620009	UDP send	0xFF800702B400016E0080
UDP send	0xFF800702B40004620009	UDP send	0xFF800702B400026E0080

UDP send	0xFF800702B400036E0080	UDP send	0xFF800702B40002C80001
UDP send	0xFF800702B400046E0080	UDP send	0xFF800702B40003C80001
UDP send	0xFF800702B400816E0080	UDP send	0xFF800702B40004C80001
UDP send	0xFF800702B400826E0080	UDP send	0xFF800702B40081C80001
UDP send	0xFF800702B400836E0080	UDP send	0xFF800702B40082C80001
UDP send	0xFF800702B400846E0080	UDP send	0xFF800702B40083C80001
UDP send	0xFF800702B40000600000	UDP send	0xFF800702B40084C80001
UDP send	0xFF800702B40001C00001	UDP send	0xFF800702B400010E0001
UDP send	0xFF800702B40002C00001	UDP send	0xFF800702B400020E0000
UDP send	0xFF800702B40003C00001	UDP send	0xFF800702B400030E0001
UDP send	0xFF800702B40004C00001	UDP send	0xFF800702B400040E0000
UDP send	0xFF800702B40081C00001	UDP send	0xFF800702B400810E0001
UDP send	0xFF800702B40082C00001	UDP send	0xFF800702B400820E0000
UDP send	0xFF800702B40083C00001	UDP send	0xFF800702B400830E0001
UDP send	0xFF800702B40084C00001	UDP send	0xFF800702B400840E0000
UDP send	0xFF800702B40001C20000	UDP send	0xFF800702B40001700800
UDP send	0xFF800702B40002C20100	UDP send	0xFF800702B40002700800
UDP send	0xFF800702B40003C20200	UDP send	0xFF800702B40003700800
UDP send	0xFF800702B40004C20300	UDP send	0xFF800702B40004700800
UDP send	0xFF800702B40005C20400	UDP send	0xFF800702B40081700800
UDP send	0xFF800702B40006C20500	UDP send	0xFF800702B40082700800
UDP send	0xFF800702B40007C20600	UDP send	0xFF800702B40083700800
UDP send	0xFF800702B40008C20700	UDP send	0xFF800702B40084700800
UDP send	0xFF800702B40081C20000	UDP send	0xFF800702B40001B00001
UDP send	0xFF800702B40082C20100	UDP send	0xFF800702B40002B00001
UDP send	0xFF800702B40083C20000	UDP send	0xFF800702B40003B00001
UDP send	0xFF800702B40084C20100	UDP send	0xFF800702B40004B00001
UDP send	0xFF800702B40085C20000	UDP send	0xFF800702B40081B00001
UDP send	0xFF800702B40086C20100	UDP send	0xFF800702B40082B00001
UDP send	0xFF800702B40087C20000	UDP send	0xFF800702B40083B00001
UDP send	0xFF800702B40088C20100	UDP send	0xFF800702B40084B00001
UDP send	0xFF800702B40001C60002	UDP send	0xFF800702B40001B400EB
UDP send	0xFF800702B40002C60002	UDP send	0xFF800702B40002B400E8
UDP send	0xFF800702B40003C60002	UDP send	0xFF800702B40003B400E4
UDP send	0xFF800702B40004C60002	UDP send	0xFF800702B40004B400DA
UDP send	0xFF800702B40081C60002	UDP send	0xFF800702B40005B400F0
UDP send	0xFF800702B40082C60002	UDP send	0xFF800702B40006B400E8
UDP send	0xFF800702B40083C60002	UDP send	0xFF800702B40007B400EB
UDP send	0xFF800702B40084C60002	UDP send	0xFF800702B40008B400E6
UDP send	0xFF800702B40001C80001	UDP send	0xFF800702B40081B40000

UDP send	0xFF800702B40082B40000	UDP send	0xFF800702B40002CE003C
UDP send	0xFF800702B40083B40000	UDP send	0xFF800702B40003CE0000
UDP send	0xFF800702B40084B40000	UDP send	0xFF800702B40004CE0032
UDP send	0xFF800702B40085B40000	UDP send	0xFF800702B40005CE0014
UDP send	0xFF800702B40086B40000	UDP send	0xFF800702B40006CE0000
UDP send	0xFF800702B40087B40000	UDP send	0xFF800702B40007CE0000
UDP send	0xFF800702B40088B40000	UDP send	0xFF800702B40008CE000A
UDP send	0xFF800702B40001B600DD	UDP send	0xFF800702B40081CE0000
UDP send	0xFF800702B40002B600DB	UDP send	0xFF800702B40082CE0000
UDP send	0xFF800702B40003B600E6	UDP send	0xFF800702B40083CE0000
UDP send	0xFF800702B40004B600E6	UDP send	0xFF800702B40084CE0000
UDP send	0xFF800702B40005B600E1	UDP send	0xFF800702B40085CE0000
UDP send	0xFF800702B40006B600E6	UDP send	0xFF800702B40086CE0000
UDP send	0xFF800702B40007B600DD	UDP send	0xFF800702B40087CE0000
UDP send	0xFF800702B40008B600F5	UDP send	0xFF800702B40088CE0000
UDP send	0xFF800702B40081B60000	UDP send	0xFF800702B40001B8000E
UDP send	0xFF800702B40082B60000	UDP send	0xFF800702B40002B8000E
UDP send	0xFF800702B40083B60000	UDP send	0xFF800702B40003B8000E
UDP send	0xFF800702B40084B60000	UDP send	0xFF800702B40004B8000E
UDP send	0xFF800702B40085B60000	UDP send	0xFF800702B40005B8000E
UDP send	0xFF800702B40086B60000	UDP send	0xFF800702B40006B8000E
UDP send	0xFF800702B40087B60000	UDP send	0xFF800702B40007B8000E
UDP send	0xFF800702B40088B60000	UDP send	0xFF800702B40008B8000E
UDP send	0xFF800702B40001CC0050	UDP send	0xFF800702B40081B80000
UDP send	0xFF800702B40002CC0064	UDP send	0xFF800702B40082B80000
UDP send	0xFF800702B40003CC000A	UDP send	0xFF800702B40083B80000
UDP send	0xFF800702B40004CC000A	UDP send	0xFF800702B40084B80000
UDP send	0xFF800702B40005CC0000	UDP send	0xFF800702B40085B80000
UDP send	0xFF800702B40006CC0000	UDP send	0xFF800702B40086B80000
UDP send	0xFF800702B40007CC0000	UDP send	0xFF800702B40087B80000
UDP send	0xFF800702B40008CC0000	UDP send	0xFF800702B40088B80000
UDP send	0xFF800702B40081CC0000	UDP send	0xFF800702B40001BA000E
UDP send	0xFF800702B40082CC0000	UDP send	0xFF800702B40002BA000E
UDP send	0xFF800702B40083CC0000	UDP send	0xFF800702B40003BA000E
UDP send	0xFF800702B40084CC0000	UDP send	0xFF800702B40004BA000E
UDP send	0xFF800702B40085CC0000	UDP send	0xFF800702B40005BA000E
UDP send	0xFF800702B40086CC0000	UDP send	0xFF800702B40006BA000E
UDP send	0xFF800702B40087CC0000	UDP send	0xFF800702B40007BA000E
UDP send	0xFF800702B40088CC0000	UDP send	0xFF800702B40008BA000E
UDP send	0xFF800702B40001CE0000	UDP send	0xFF800702B40081BA0000

UDP send 0xFF800702B40082BA0000
UDP send 0xFF800702B40083BA0000
UDP send 0xFF800702B40084BA0000
UDP send 0xFF800702B40085BA0000
UDP send 0xFF800702B40086BA0000
UDP send 0xFF800702B40087BA0000
UDP send 0xFF800702B40088BA0000
UDP send 0xFF800702B40001BC007B
UDP send 0xFF800702B40002BC007E
UDP send 0xFF800702B40003BC0088
UDP send 0xFF800702B40004BC0080
UDP send 0xFF800702B40005BC0080
UDP send 0xFF800702B40006BC008E
UDP send 0xFF800702B40007BC0085
UDP send 0xFF800702B40008BC007D
UDP send 0xFF800702B40081BC0000
UDP send 0xFF800702B40082BC0000
UDP send 0xFF800702B40083BC0000
UDP send 0xFF800702B40084BC0000
UDP send 0xFF800702B40085BC0000
UDP send 0xFF800702B40086BC0000
UDP send 0xFF800702B40087BC0000
UDP send 0xFF800702B40088BC0000
UDP send 0xFF800702B40001BE0085
UDP send 0xFF800702B40002BE0082
UDP send 0xFF800702B40003BE0078
UDP send 0xFF800702B40004BE0080
UDP send 0xFF800702B40005BE0080
UDP send 0xFF800702B40006BE0072
UDP send 0xFF800702B40007BE007B
UDP send 0xFF800702B40008BE0083
UDP send 0xFF800702B40081BE0000
UDP send 0xFF800702B40082BE0000
UDP send 0xFF800702B40083BE0000
UDP send 0xFF800702B40084BE0000
UDP send 0xFF800702B40085BE0000
UDP send 0xFF800702B40086BE0000
UDP send 0xFF800702B40087BE0000
UDP send 0xFF800702B40088BE0000
UDP send 0xFF800702B40001D80005

UDP send 0xFF800702B40002D80005
UDP send 0xFF800702B40003D80005
UDP send 0xFF800702B40004D80005
UDP send 0xFF800702B40081D80005
UDP send 0xFF800702B40082D80005
UDP send 0xFF800702B40083D80005
UDP send 0xFF800702B40084D80005
UDP send 0xFF800702B40001DA0005
UDP send 0xFF800702B40002DA0005
UDP send 0xFF800702B40003DA0005
UDP send 0xFF800702B40004DA0005
UDP send 0xFF800702B40081DA0005
UDP send 0xFF800702B40082DA0005
UDP send 0xFF800702B40083DA0005
UDP send 0xFF800702B40084DA0005
UDP send 0xFF800702B40001DC0017
UDP send 0xFF800702B40002DC0017
UDP send 0xFF800702B40003DC0017
UDP send 0xFF800702B40004DC0017
UDP send 0xFF800702B40081DC0017
UDP send 0xFF800702B40082DC0017
UDP send 0xFF800702B40083DC0017
UDP send 0xFF800702B40084DC0017
UDP send 0xFF800702B40001DE0000
UDP send 0xFF800702B40002DE0000
UDP send 0xFF800702B40003DE0000
UDP send 0xFF800702B40004DE0000
UDP send 0xFF800702B40081DE0000
UDP send 0xFF800702B40082DE0000
UDP send 0xFF800702B40083DE0000
UDP send 0xFF800702B40084DE0000
UDP send 0xFF800702B40001100032
UDP send 0xFF800702B40002100032
UDP send 0xFF800702B40003100032
UDP send 0xFF800702B40004100032
UDP send 0xFF800702B40081100032
UDP send 0xFF800702B40082100032
UDP send 0xFF800702B40083100032
UDP send 0xFF800702B40084100032
UDP send 0xFF800702B40001D00000

UDP send 0xFF800702B40002D00000
UDP send 0xFF800702B40003D00000
UDP send 0xFF800702B40004D00000
UDP send 0xFF800702B40081D00000
UDP send 0xFF800702B40082D00000
UDP send 0xFF800702B40083D00000
UDP send 0xFF800702B40084D00000
UDP send 0xFF800702B40040900000
UDP send 0xFF800702B40040900001
UDP send 0xFF800702B40040900000
UDP send 0xFF800702B40040020000
UDP send 0xFF800702B4000174000A
UDP send 0xFF800702B4000274000A
UDP send 0xFF800702B4000374000A
UDP send 0xFF800702B4000474000A
UDP send 0xFF800702B4008174000A
UDP send 0xFF800702B4008274000A
UDP send 0xFF800702B4008374000A
UDP send 0xFF800702B4008474000A
UDP send 0xFF800702B4000172000F
UDP send 0xFF800702B4000272000F
UDP send 0xFF800702B4000372000F
UDP send 0xFF800702B4000472000F
UDP send 0xFF800702B4008172000F
UDP send 0xFF800702B4008272000F
UDP send 0xFF800702B4008372000F
UDP send 0xFF800702B4008472000F
UDP send 0xFF800702B40040260000
UDP send 0xFF800702B40001760000
UDP send 0xFF800702B40002760000
UDP send 0xFF800702B40003760000
UDP send 0xFF800702B40004760000
UDP send 0xFF800702B40081760000
UDP send 0xFF800702B40082760000
UDP send 0xFF800702B40083760000
UDP send 0xFF800702B40084760000
UDP send 0xFF800702B40001780000
UDP send 0xFF800702B40002780000
UDP send 0xFF800702B40003780000
UDP send 0xFF800702B40004780000
UDP send 0xFF800702B40081780000
UDP send 0xFF800702B40082780000
UDP send 0xFF800702B40083780000
UDP send 0xFF800702B40084780000

UDP send 0xFF800702B40001EA0014
 UDP send 0xFF800702B40002EA0014
 UDP send 0xFF800702B40003EA0014
 UDP send 0xFF800702B40004EA0014
 UDP send 0xFF800702B40081EA0014
 UDP send 0xFF800702B40082EA0014
 UDP send 0xFF800702B40083EA0014
 UDP send 0xFF800702B40084EA0014
 UDP send 0xFF800702B40001EC000A
 UDP send 0xFF800702B40002EC000A
 UDP send 0xFF800702B40003EC000A
 UDP send 0xFF800702B40004EC000A
 UDP send 0xFF800702B40081EC000A
 UDP send 0xFF800702B40082EC000A
 UDP send 0xFF800702B40083EC000A
 UDP send 0xFF800702B40084EC000A
 UDP send 0xFF800702B40001EE0014
 UDP send 0xFF800702B40002EE0014
 UDP send 0xFF800702B40003EE0014
 UDP send 0xFF800702B40004EE0014
 UDP send 0xFF800702B40081EE0014
 UDP send 0xFF800702B40082EE0014
 UDP send 0xFF800702B40083EE0014
 UDP send 0xFF800702B40084EE0014
 UDP send 0xFF800702B40001D60000
 UDP send 0xFF800702B40002D60000
 UDP send 0xFF800702B40003D60000
 UDP send 0xFF800702B40004D60000
 UDP send 0xFF800702B40081D60000
 UDP send 0xFF800702B40082D60000
 UDP send 0xFF800702B40083D60000
 UDP send 0xFF800702B40084D60000
 UDP send 0xFF800702B40001820001
 UDP send 0xFF800702B40002820001
 UDP send 0xFF800702B40003820001
 UDP send 0xFF800702B40004820001
 UDP send 0xFF800702B40081820001
 UDP send 0xFF800702B40082820001
 UDP send 0xFF800702B40083820001
 UDP send 0xFF800702B40084820001

UDP send 0xFF800702B400017C0001
 UDP send 0xFF800702B400027C0001
 UDP send 0xFF800702B400037C0001
 UDP send 0xFF800702B400047C0001
 UDP send 0xFF800702B400817C0001
 UDP send 0xFF800702B400827C0001
 UDP send 0xFF800702B400837C0001
 UDP send 0xFF800702B400847C0001
 UDP send 0xFF800702B40001860005
 UDP send 0xFF800702B40002860005
 UDP send 0xFF800702B40003860005
 UDP send 0xFF800702B40004860005
 UDP send 0xFF800702B40081860005
 UDP send 0xFF800702B40082860005
 UDP send 0xFF800702B40083860005
 UDP send 0xFF800702B40084860005
 UDP send 0xFF800702B40001880004
 UDP send 0xFF800702B40002880004
 UDP send 0xFF800702B40003880005
 UDP send 0xFF800702B40004880005
 UDP send 0xFF800702B40081880005
 UDP send 0xFF800702B40082880005
 UDP send 0xFF800702B40083880005
 UDP send 0xFF800702B40084880005
 UDP send 0xFF800702B40040720000
 UDP send 0xFF800702B40001400000
 UDP send 0xFF800702B40001400001
 UDP send 0xFF800702B40001400000
 UDP send 0xFF800702B40081400000
 UDP send 0xFF800702B40081400001
 UDP send 0xFF800702B40081400000

2) Config setting command list example

UDP send	0xFF800702B4004000001	UDP send	0xFF800702B40081620009
UDP send	0xFF800702B40040060000	UDP send	0xFF800702B40082620009
UDP send	0xFF800702B40040080000	UDP send	0xFF800702B40083620009
UDP send	0xFF800702B400400A2540	UDP send	0xFF800702B40084620009
UDP send	0xFF800702B400400CBE40	UDP send	0xFF800702B40001640019
UDP send	0xFF800702B400402E0100	UDP send	0xFF800702B40002640019
UDP send	0xFF800702B40040300001	UDP send	0xFF800702B40003640019
UDP send	0xFF800702B400403200FF	UDP send	0xFF800702B40004640019
UDP send	0xFF800702B400008C0007	UDP send	0xFF800702B40081640019
UDP send	0xFF800702B400011A0001	UDP send	0xFF800702B40082640019
UDP send	0xFF800702B400021A0001	UDP send	0xFF800702B40083640019
UDP send	0xFF800702B400031A0001	UDP send	0xFF800702B40084640019
UDP send	0xFF800702B400041A0001	UDP send	0xFF800702B4000166001E
UDP send	0xFF800702B400811A0001	UDP send	0xFF800702B4000266001E
UDP send	0xFF800702B400821A0001	UDP send	0xFF800702B4000366001E
UDP send	0xFF800702B400831A0001	UDP send	0xFF800702B4000466001E
UDP send	0xFF800702B400841A0001	UDP send	0xFF800702B4008166001E
UDP send	0xFF800702B40040360000	UDP send	0xFF800702B4008266001E
UDP send	0xFF800702B400010C0004	UDP send	0xFF800702B4008366001E
UDP send	0xFF800702B400020C0004	UDP send	0xFF800702B4008466001E
UDP send	0xFF800702B400030C0004	UDP send	0xFF800702B4000168001E
UDP send	0xFF800702B400040C0004	UDP send	0xFF800702B4000268001E
UDP send	0xFF800702B400810C0004	UDP send	0xFF800702B4000368001E
UDP send	0xFF800702B400820C0004	UDP send	0xFF800702B4000468001E
UDP send	0xFF800702B400830C0004	UDP send	0xFF800702B4008168001E
UDP send	0xFF800702B400840C0004	UDP send	0xFF800702B4008268001E
UDP send	0xFF800702B40001600007	UDP send	0xFF800702B4008368001E
UDP send	0xFF800702B40002600007	UDP send	0xFF800702B4008468001E
UDP send	0xFF800702B40003600007	UDP send	0xFF800702B400016A1F40
UDP send	0xFF800702B40004600007	UDP send	0xFF800702B400026A1F40
UDP send	0xFF800702B40081600007	UDP send	0xFF800702B400036A1F40
UDP send	0xFF800702B40082600007	UDP send	0xFF800702B400046A1F40
UDP send	0xFF800702B40083600007	UDP send	0xFF800702B400816A1F40
UDP send	0xFF800702B40084600007	UDP send	0xFF800702B400826A1F40
UDP send	0xFF800702B40001620009	UDP send	0xFF800702B400836A1F40
UDP send	0xFF800702B40002620009	UDP send	0xFF800702B400846A1F40
UDP send	0xFF800702B40003620009	UDP send	0xFF800702B400016E0080
UDP send	0xFF800702B40004620009	UDP send	0xFF800702B400026E0080

UDP send	0xFF800702B400036E0080	UDP send	0xFF800702B40002C80001
UDP send	0xFF800702B400046E0080	UDP send	0xFF800702B40003C80001
UDP send	0xFF800702B400816E0080	UDP send	0xFF800702B40004C80001
UDP send	0xFF800702B400826E0080	UDP send	0xFF800702B40081C80001
UDP send	0xFF800702B400836E0080	UDP send	0xFF800702B40082C80001
UDP send	0xFF800702B400846E0080	UDP send	0xFF800702B40083C80001
UDP send	0xFF800702B40000600000	UDP send	0xFF800702B40084C80001
UDP send	0xFF800702B40001C00001	UDP send	0xFF800702B400010E0001
UDP send	0xFF800702B40002C00001	UDP send	0xFF800702B400020E0000
UDP send	0xFF800702B40003C00001	UDP send	0xFF800702B400030E0001
UDP send	0xFF800702B40004C00001	UDP send	0xFF800702B400040E0000
UDP send	0xFF800702B40081C00001	UDP send	0xFF800702B400810E0001
UDP send	0xFF800702B40082C00001	UDP send	0xFF800702B400820E0000
UDP send	0xFF800702B40083C00001	UDP send	0xFF800702B400830E0001
UDP send	0xFF800702B40084C00001	UDP send	0xFF800702B400840E0000
UDP send	0xFF800702B40001C20000	UDP send	0xFF800702B40001700800
UDP send	0xFF800702B40002C20100	UDP send	0xFF800702B40002700800
UDP send	0xFF800702B40003C20200	UDP send	0xFF800702B40003700800
UDP send	0xFF800702B40004C20300	UDP send	0xFF800702B40004700800
UDP send	0xFF800702B40005C20400	UDP send	0xFF800702B40081700800
UDP send	0xFF800702B40006C20500	UDP send	0xFF800702B40082700800
UDP send	0xFF800702B40007C20600	UDP send	0xFF800702B40083700800
UDP send	0xFF800702B40008C20700	UDP send	0xFF800702B40084700800
UDP send	0xFF800702B40081C20000	UDP send	0xFF800702B40001B00001
UDP send	0xFF800702B40082C20100	UDP send	0xFF800702B40002B00001
UDP send	0xFF800702B40083C20000	UDP send	0xFF800702B40003B00001
UDP send	0xFF800702B40084C20100	UDP send	0xFF800702B40004B00001
UDP send	0xFF800702B40085C20000	UDP send	0xFF800702B40081B00001
UDP send	0xFF800702B40086C20100	UDP send	0xFF800702B40082B00001
UDP send	0xFF800702B40087C20000	UDP send	0xFF800702B40083B00001
UDP send	0xFF800702B40088C20100	UDP send	0xFF800702B40084B00001
UDP send	0xFF800702B40001C60002	UDP send	0xFF800702B40001B400E8
UDP send	0xFF800702B40002C60002	UDP send	0xFF800702B40002B400E8
UDP send	0xFF800702B40003C60002	UDP send	0xFF800702B40003B400E4
UDP send	0xFF800702B40004C60002	UDP send	0xFF800702B40004B400DA
UDP send	0xFF800702B40081C60002	UDP send	0xFF800702B40005B400F0
UDP send	0xFF800702B40082C60002	UDP send	0xFF800702B40006B400E8
UDP send	0xFF800702B40083C60002	UDP send	0xFF800702B40007B400EB
UDP send	0xFF800702B40084C60002	UDP send	0xFF800702B40008B400E6
UDP send	0xFF800702B40001C80001	UDP send	0xFF800702B40081B40000

UDP send	0xFF800702B40082B40000	UDP send	0xFF800702B40002CE003C
UDP send	0xFF800702B40083B40000	UDP send	0xFF800702B40003CE0000
UDP send	0xFF800702B40084B40000	UDP send	0xFF800702B40004CE0032
UDP send	0xFF800702B40085B40000	UDP send	0xFF800702B40005CE0014
UDP send	0xFF800702B40086B40000	UDP send	0xFF800702B40006CE0000
UDP send	0xFF800702B40087B40000	UDP send	0xFF800702B40007CE0000
UDP send	0xFF800702B40088B40000	UDP send	0xFF800702B40008CE000A
UDP send	0xFF800702B40001B600DD	UDP send	0xFF800702B40081CE0000
UDP send	0xFF800702B40002B600DB	UDP send	0xFF800702B40082CE0000
UDP send	0xFF800702B40003B600E6	UDP send	0xFF800702B40083CE0000
UDP send	0xFF800702B40004B600E6	UDP send	0xFF800702B40084CE0000
UDP send	0xFF800702B40005B600E1	UDP send	0xFF800702B40085CE0000
UDP send	0xFF800702B40006B600E6	UDP send	0xFF800702B40086CE0000
UDP send	0xFF800702B40007B600DD	UDP send	0xFF800702B40087CE0000
UDP send	0xFF800702B40008B600F5	UDP send	0xFF800702B40088CE0000
UDP send	0xFF800702B40081B60000	UDP send	0xFF800702B40001B8000E
UDP send	0xFF800702B40082B60000	UDP send	0xFF800702B40002B8000E
UDP send	0xFF800702B40083B60000	UDP send	0xFF800702B40003B8000E
UDP send	0xFF800702B40084B60000	UDP send	0xFF800702B40004B8000E
UDP send	0xFF800702B40085B60000	UDP send	0xFF800702B40005B8000E
UDP send	0xFF800702B40086B60000	UDP send	0xFF800702B40006B8000E
UDP send	0xFF800702B40087B60000	UDP send	0xFF800702B40007B8000E
UDP send	0xFF800702B40088B60000	UDP send	0xFF800702B40008B8000E
UDP send	0xFF800702B40001CC0050	UDP send	0xFF800702B40081B80000
UDP send	0xFF800702B40002CC0064	UDP send	0xFF800702B40082B80000
UDP send	0xFF800702B40003CC000A	UDP send	0xFF800702B40083B80000
UDP send	0xFF800702B40004CC000A	UDP send	0xFF800702B40084B80000
UDP send	0xFF800702B40005CC0000	UDP send	0xFF800702B40085B80000
UDP send	0xFF800702B40006CC0000	UDP send	0xFF800702B40086B80000
UDP send	0xFF800702B40007CC0000	UDP send	0xFF800702B40087B80000
UDP send	0xFF800702B40008CC0000	UDP send	0xFF800702B40088B80000
UDP send	0xFF800702B40081CC0000	UDP send	0xFF800702B40001BA000E
UDP send	0xFF800702B40082CC0000	UDP send	0xFF800702B40002BA000E
UDP send	0xFF800702B40083CC0000	UDP send	0xFF800702B40003BA000E
UDP send	0xFF800702B40084CC0000	UDP send	0xFF800702B40004BA000E
UDP send	0xFF800702B40085CC0000	UDP send	0xFF800702B40005BA000E
UDP send	0xFF800702B40086CC0000	UDP send	0xFF800702B40006BA000E
UDP send	0xFF800702B40087CC0000	UDP send	0xFF800702B40007BA000E
UDP send	0xFF800702B40088CC0000	UDP send	0xFF800702B40008BA000E
UDP send	0xFF800702B40001CE0000	UDP send	0xFF800702B40081BA0000

UDP send	0xFF800702B40082BA0000	UDP send	0xFF800702B40002D80005
UDP send	0xFF800702B40083BA0000	UDP send	0xFF800702B40003D80005
UDP send	0xFF800702B40084BA0000	UDP send	0xFF800702B40004D80005
UDP send	0xFF800702B40085BA0000	UDP send	0xFF800702B40081D80005
UDP send	0xFF800702B40086BA0000	UDP send	0xFF800702B40082D80005
UDP send	0xFF800702B40087BA0000	UDP send	0xFF800702B40083D80005
UDP send	0xFF800702B40088BA0000	UDP send	0xFF800702B40084D80005
UDP send	0xFF800702B40001BC007B	UDP send	0xFF800702B40001DA0005
UDP send	0xFF800702B40002BC007E	UDP send	0xFF800702B40002DA0005
UDP send	0xFF800702B40003BC0088	UDP send	0xFF800702B40003DA0005
UDP send	0xFF800702B40004BC0080	UDP send	0xFF800702B40004DA0005
UDP send	0xFF800702B40005BC0080	UDP send	0xFF800702B40081DA0005
UDP send	0xFF800702B40006BC008E	UDP send	0xFF800702B40082DA0005
UDP send	0xFF800702B40007BC0085	UDP send	0xFF800702B40083DA0005
UDP send	0xFF800702B40008BC007D	UDP send	0xFF800702B40084DA0005
UDP send	0xFF800702B40081BC0000	UDP send	0xFF800702B40001DC0017
UDP send	0xFF800702B40082BC0000	UDP send	0xFF800702B40002DC0017
UDP send	0xFF800702B40083BC0000	UDP send	0xFF800702B40003DC0017
UDP send	0xFF800702B40084BC0000	UDP send	0xFF800702B40004DC0017
UDP send	0xFF800702B40085BC0000	UDP send	0xFF800702B40081DC0017
UDP send	0xFF800702B40086BC0000	UDP send	0xFF800702B40082DC0017
UDP send	0xFF800702B40087BC0000	UDP send	0xFF800702B40083DC0017
UDP send	0xFF800702B40088BC0000	UDP send	0xFF800702B40084DC0017
UDP send	0xFF800702B40001BE0085	UDP send	0xFF800702B40001DE0000
UDP send	0xFF800702B40002BE0082	UDP send	0xFF800702B40002DE0000
UDP send	0xFF800702B40003BE0078	UDP send	0xFF800702B40003DE0000
UDP send	0xFF800702B40004BE0080	UDP send	0xFF800702B40004DE0000
UDP send	0xFF800702B40005BE0080	UDP send	0xFF800702B40081DE0000
UDP send	0xFF800702B40006BE0072	UDP send	0xFF800702B40082DE0000
UDP send	0xFF800702B40007BE007B	UDP send	0xFF800702B40083DE0000
UDP send	0xFF800702B40008BE0083	UDP send	0xFF800702B40084DE0000
UDP send	0xFF800702B40081BE0000	UDP send	0xFF800702B40001100032
UDP send	0xFF800702B40082BE0000	UDP send	0xFF800702B40002100032
UDP send	0xFF800702B40083BE0000	UDP send	0xFF800702B40003100032
UDP send	0xFF800702B40084BE0000	UDP send	0xFF800702B40004100032
UDP send	0xFF800702B40085BE0000	UDP send	0xFF800702B40081100032
UDP send	0xFF800702B40086BE0000	UDP send	0xFF800702B40082100032
UDP send	0xFF800702B40087BE0000	UDP send	0xFF800702B40083100032
UDP send	0xFF800702B40088BE0000	UDP send	0xFF800702B40084100032
UDP send	0xFF800702B40001D80005	UDP send	0xFF800702B40001D00000

UDP send	0xFF800702B40002D00000	UDP send	0xFF800702B40081780000
UDP send	0xFF800702B40003D00000	UDP send	0xFF800702B40082780000
UDP send	0xFF800702B40004D00000	UDP send	0xFF800702B40083780000
UDP send	0xFF800702B40081D00000	UDP send	0xFF800702B40084780000
UDP send	0xFF800702B40082D00000	UDP send	0xFF800702B400017A0010
UDP send	0xFF800702B40083D00000	UDP send	0xFF800702B400027A0010
UDP send	0xFF800702B40084D00000	UDP send	0xFF800702B400037A0010
UDP send	0xFF800702B40040900000	UDP send	0xFF800702B400047A0010
UDP send	0xFF800702B40040900001	UDP send	0xFF800702B400817A0010
UDP send	0xFF800702B40040900000	UDP send	0xFF800702B400827A0010
UDP send	0xFF800702B40040020000	UDP send	0xFF800702B400837A0010
UDP send	0xFF800702B4000174000A	UDP send	0xFF800702B400847A0010
UDP send	0xFF800702B4000274000A	UDP send	0xFF800702B40040700019
UDP send	0xFF800702B4000374000A	UDP send	0xFF800702B40001800000
UDP send	0xFF800702B4000474000A	UDP send	0xFF800702B40002800000
UDP send	0xFF800702B4008174000A	UDP send	0xFF800702B40003800000
UDP send	0xFF800702B4008274000A	UDP send	0xFF800702B40004800000
UDP send	0xFF800702B4008374000A	UDP send	0xFF800702B40081800000
UDP send	0xFF800702B4008474000A	UDP send	0xFF800702B40082800000
UDP send	0xFF800702B4000172000F	UDP send	0xFF800702B40083800000
UDP send	0xFF800702B4000272000F	UDP send	0xFF800702B40084800000
UDP send	0xFF800702B4000372000F	UDP send	0xFF800702B40000480001
UDP send	0xFF800702B4000472000F	UDP send	0xFF800702B400404A0008
UDP send	0xFF800702B4008172000F	UDP send	0xFF800702B40001840037
UDP send	0xFF800702B4008272000F	UDP send	0xFF800702B40002840023
UDP send	0xFF800702B4008372000F	UDP send	0xFF800702B4000384003C
UDP send	0xFF800702B4008472000F	UDP send	0xFF800702B4000484003E
UDP send	0xFF800702B40040260000	UDP send	0xFF800702B4008184004B
UDP send	0xFF800702B40001760000	UDP send	0xFF800702B4008284004E
UDP send	0xFF800702B40002760000	UDP send	0xFF800702B40083840046
UDP send	0xFF800702B40003760000	UDP send	0xFF800702B40084840049
UDP send	0xFF800702B40004760000	UDP send	0xFF800702B400008A0000
UDP send	0xFF800702B40081760000	UDP send	0xFF800702B40001E8000A
UDP send	0xFF800702B40082760000	UDP send	0xFF800702B40002E8000A
UDP send	0xFF800702B40083760000	UDP send	0xFF800702B40003E8000A
UDP send	0xFF800702B40084760000	UDP send	0xFF800702B40004E8000A
UDP send	0xFF800702B40001780000	UDP send	0xFF800702B40081E8000A
UDP send	0xFF800702B40002780000	UDP send	0xFF800702B40082E8000A
UDP send	0xFF800702B40003780000	UDP send	0xFF800702B40083E8000A
UDP send	0xFF800702B40004780000	UDP send	0xFF800702B40084E8000A

UDP send	0xFF800702B40001EA0014	UDP send	0xFF800702B40003820001
UDP send	0xFF800702B40002EA0014	UDP send	0xFF800702B40004820001
UDP send	0xFF800702B40003EA0014	UDP send	0xFF800702B40081820001
UDP send	0xFF800702B40004EA0014	UDP send	0xFF800702B40082820001
UDP send	0xFF800702B40081EA0014	UDP send	0xFF800702B40083820001
UDP send	0xFF800702B40082EA0014	UDP send	0xFF800702B40084820001
UDP send	0xFF800702B40083EA0014	UDP send	0xFF800702B400017C0001
UDP send	0xFF800702B40084EA0014	UDP send	0xFF800702B400027C0001
UDP send	0xFF800702B40001EC000A	UDP send	0xFF800702B400037C0001
UDP send	0xFF800702B40002EC000A	UDP send	0xFF800702B400047C0001
UDP send	0xFF800702B40003EC000A	UDP send	0xFF800702B400817C0001
UDP send	0xFF800702B40004EC000A	UDP send	0xFF800702B400827C0001
UDP send	0xFF800702B40081EC000A	UDP send	0xFF800702B400837C0001
UDP send	0xFF800702B40082EC000A	UDP send	0xFF800702B400847C0001
UDP send	0xFF800702B40083EC000A	UDP send	0xFF800702B40001860005
UDP send	0xFF800702B40084EC000A	UDP send	0xFF800702B40002860005
UDP send	0xFF800702B40001EE0014	UDP send	0xFF800702B40003860005
UDP send	0xFF800702B40002EE0014	UDP send	0xFF800702B40004860005
UDP send	0xFF800702B40003EE0014	UDP send	0xFF800702B40081860005
UDP send	0xFF800702B40004EE0014	UDP send	0xFF800702B40082860005
UDP send	0xFF800702B40081EE0014	UDP send	0xFF800702B40083860005
UDP send	0xFF800702B40082EE0014	UDP send	0xFF800702B40084860005
UDP send	0xFF800702B40083EE0014	UDP send	0xFF800702B40001880004
UDP send	0xFF800702B40084EE0014	UDP send	0xFF800702B40002880004
UDP send	0xFF800702B40001D60000	UDP send	0xFF800702B40003880005
UDP send	0xFF800702B40002D60000	UDP send	0xFF800702B40004880005
UDP send	0xFF800702B40003D60000	UDP send	0xFF800702B40081880005
UDP send	0xFF800702B40004D60000	UDP send	0xFF800702B40082880005
UDP send	0xFF800702B40081D60000	UDP send	0xFF800702B40083880005
UDP send	0xFF800702B40082D60000	UDP send	0xFF800702B40084880005
UDP send	0xFF800702B40083D60000	UDP send	0xFF800702B40040720000
UDP send	0xFF800702B40084D60000		
UDP send	0xFF800702B40001820001		
UDP send	0xFF800702B40002820001		

9. End

Click "quit" in "File" of the menu bar. After clicking, this application ends, and the screen disappears. The next time you start up, the settings at the end will be reflected.

CONTACT INFORMATION

TechnoAP Co., Ltd.	
ADD	2976-15 Mawatari, Hitachinaka-shi, Ibaraki, 312-0012, Japan
TEL	+81-29-350-8011
FAX	+81-29-352-9013
WEB	http://www.techno-ap.com
Email	order@techno-ap.com

[Business Hour]: Monday to Friday, 9:30 to 17:00 (Japan Time)